

rea | *Azalea Max*



A tailor-made solution for heavier users

The *Rea Azalea Max* is designed to meet the specific needs of disabled users with extra weight. As a specific good option we recommend power tilt and backrest recline. The *Rea Azalea Max* is a niche product and a complementary to the existing *Azalea* family from Invacare. Now we could offer *Azalea Max* with additional frame colours.





Features and Options

Continued Compatibility

The *Rea Azalea Max* offers the same features as the rest of the *Rea Azalea* family such as tilt-in-space with weight-shifting and DSS (Dual Stability System) mechanism. Stability in this wheelchair originates from body length compensation in the back and leg rests, beneficial for users with limited movement and that have difficulty changing position.

High ergonomic push handle

The ergonomically shaped push handles are positioned to optimize the pushing for the attendant. An alternate version with higher handles offers a different position for the hands and is available as an option.



Two-way backrest

The lower and upper parts of the shaped backrest are reversible, giving users the option of more space in the seat or chest area.



Footplate converter

Divided into three parts, the footplate comes with a removable mid section, to facilitate transfers.



Power tilt and recline

With the hand control the user or the attendant can easily change seating position.

Technical data

Azalea Max	550; 610 + 50 mm	530 - 600 mm	400 / 450 mm	500 + 80 mm	320 - 420 mm	330 - 500 mm	1° - 20°
Azalea Max Assist	550; 610 + 100 mm	530 - 600 mm	400 / 450 mm	500 + 80 mm	320 - 420 mm	330 - 500 mm	1° - 20°
Azalea Max	0°-30°	SW + 250 mm	1020 - 1240 mm	1050 - 1120 mm	54 kg		
Azalea Max Assist	0°-30°	SW + 220 mm	1020 - 1240 mm	1050 - 1120 mm	54 kg		
		For other <i>Rea Azalea</i> accessories, please consult the <i>Rea Azalea</i> brochure.					
Azalea Max	33 kg	The <i>Rea Azalea Max</i> complies with ISO 7176-19 specifications. Note that special conditions apply regarding accessories. For detailed information regarding vehicular transportation and for more pre-sale information for this wheelchair, please consult the following pages on our website: www.invacare-rea.com under Rea wheelchairs/Technical downloads.					
Azalea Max Assist	37 kg	max. 180 kg					

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© 2012 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Rea Azalea Max - EU - 02/2012




Yes, you can.®

Frame colour



Upholstery



Rea, Rea design and DSS (Dual Stability System) design are registered trademarks of Invacare International.

The *Rea Azalea Max* is designed with DSS (Dual Stability system). It provides the user with enhanced stability and safety when being tilted or the backrest is reclined.

